




#### RESEARCH ARTICLE

## Differentiating Positional Isomers of Nucleoside **Modifications by Higher-Energy Collisional Dissociation** Mass Spectrometry (HCD MS)

Manasses Jora, Andrew P. Burns, Robert L. Ross, Peter A. Lobue, Ruoxia Zhao, Cody M. Palumbo,<sup>2</sup> Peter A. Beal,<sup>2</sup> Balasubrahmanyam Addepalli,<sup>1</sup> Patrick A. Limbach<sup>1</sup>

<sup>1</sup>Rieveschl Laboratories for Mass Spectrometry, Department of Chemistry, University of Cincinnati, Cincinnati, OH 45221, USA <sup>2</sup>Department of Chemistry, University of California, Davis, California, 95616, USA



#### Ribonucleoside Positional Isomers

Abstract. The analytical identification of positional isomers (e.g., 3-,  $N^4$ -, 5-methylcytidine) within the > 160 different post-transcriptional modifications found in RNA can be challenging. Conventional liquid chromatography coupled with tandem mass spectrometry (LC-MS/MS) approaches rely on chromatographic separation for accurate identification because the collisioninduced dissociation (CID) mass spectra of these isomers nearly exclusively yield identical

analytical platform for both identification [1-5] and quantifi-

methods for modified nucleosides [11], collision-induced dis-

sociation (CID) of the protonated molecular ion (MH<sup>+</sup>) more

often than not leads to fragmentation at the N-glycosidic bond

of ribonucleosides generating a diagnostic nucleobase ion

this feature to create powerful analytical platforms for the

As noted extremely early in the development of LC-MS

nucleobase ions (BH<sub>2</sub><sup>+</sup>) from the same molecular ion (MH<sup>+</sup>). Here, we have explored higher-energy collisional dissociation (HCD) as an alternative fragmentation technique to generate more informative product ions that can be used to differentiate positional isomers. LC-MS/MS of modified nucleosides characterized using HCD led to the creation of structure- and HCD energy-specific fragmentation patterns that generated unique fingerprints, which can be used to identify individual positional isomers even when they cannot be separated chromatographically. While particularly useful for identifying positional isomers, the fingerprinting capabilities enabled by HCD also offer the potential to generate HPLC-independent spectral libraries for the rapid analysis of modified ribonucleosides.

Keywords: RNA modification, Positional isomers, Nucleoside analysis, LC-MS/MS, HCD fragmentation

Received: 17 April 2018/Revised: 22 May 2018/Accepted: 24 May 2018/Published Online: 12 June 2018

#### Introduction

ver 160 types of post-transcriptional nucleoside modifications have been documented in various types of ribonucleic acid (RNA) [1]. These chemical modifications range from simple methylation or thiolation to complex hypermodifications. They are added mostly to the nucleobase (except for ribose methylation) through simple or complex enzymatic pathways. Liquid chromatography- tandem mass spectrometry (LC-MS/MS) is now viewed as the gold standard

Electronic supplementary material The online version of this article (https:// doi.org/10.1007/s13361-018-1999-6) contains supplementary material, which is available to authorized users.

neutral loss of 146 Da). Various laboratories have exploited

identification of known and detection of previously unknown modifications [3, 8, 12–14].

cation [6–10] of modified ribonucleosides.

Correspondence to: Patrick Limbach; e-mail: limbacpa@ucmail.uc.edu

Despite the demonstrated utility of CID-based LC-MS/MS of modified nucleosides from RNA, it is well understood that chromatographic reproducibility and resolving power are necessary to deal with modified nucleosides that exist as positional isomers, as the common MH<sup>+</sup> to BH<sub>2</sub><sup>+</sup> transitions for such isomers yield equivalent mass spectral responses. Of all the ribonucleoside modifications, methylations are the most abundant class of positional isomers [1]. For example, the isomers of methylcytidine (3-,  $N^4$ -, and 5-methylcytidine) are indistinguishable based on the evaluation of m/z of MH<sup>+</sup> and BH<sub>2</sub><sup>+</sup> ions alone from the mass spectra. Such positional isomers of methylations observed in vivo have significant biological relevance. For instance, m<sup>5</sup>C, commonly found in the anticodon loop of archaeal and eukaryotic tRNAs, is linked to structural and metabolic stabilization of the tRNA [15]. Therefore, distinguishing positional isomers helps delineate the functional impact of a specific modification.

Typically, the analyst will rely on previously reported chromatographic behavior [1, 6, 11] or the use of analytical standards to identify and differentiate modified nucleosides that can be present as positional isomers. Another approach used in distinguishing positional isomers is multistage fragmentation (MS<sup>n</sup>) [16, 17], where analysis of the data generated by second and third (or more) stages of dissociation can discriminate among different structural possibilities. However, acquisition of sequential fragmentation events impacts the duty cycle of analysis, which can hinder accurate qualitative and quantitative measurements for co-eluting nucleosides. Another technique to address the challenge of ribonucleoside positional isomers is the inclusion of ion mobility separation prior to MS/MS [17]. Here, the cross-sectional area of the analyte allows for a gasphase separation, although it remains to be seen if all ribonucleoside positional isomers are amenable to separation in this manner.

Orbitrap-based mass analyzers are equipped with an ion routing multipole (IRM), a multipurpose ion storage device that has an increased peak-to-peak radiofrequency to retain a wider range of m/z values compared to traditional linear ion trap devices [18, 19]. The IRM enables an optional fragmentation mode, higher-energy collisional dissociation (HCD), during accumulation of ions [19]. Due to the increased trapping capability of the IRM, HCD MS/MS spectra yield beam-type instrument characteristics (i.e., pseudo-MS<sup>3</sup>) [20]. Such pseudo-MS<sup>3</sup> approaches provide improved efficiency for structural elucidation of small molecules [21].

Given the scientific interest in modified ribonucleosides along with the ever-increasing sample types being studied [22–24], we have been interested in establishing analytical platforms that can be used for the identification and characterization of known and unknown modifications that would be independent of (or significantly less dependent on) separation techniques and protocols. Based on the prior findings of Jensen and co-workers regarding the utility of MS<sup>n</sup> fragmentation for nucleoside identification [16] combined with the proven capabilities of HCD, in this work, we investigated whether the conventional LC-MS/MS approach for nucleoside analysis

could be enhanced by using HCD in place of CID. We found conditions for HCD that mimicked CID-type fragmentation, thereby allowing previous neutral loss techniques to still be applied during LC-MS/MS analyses. More significantly, we report that at higher energy values, HCD generates pseudo-MS<sup>n</sup> product ion spectra that can be used to differentiate positional isomers, even when those isomers would co-elute during LC-MS/MS. Building on that finding, we propose that HCD can be used to generate nucleoside-specific fingerprints, which can enable spectral library matching for any modified ribonucleoside. By using only mass spectrometry data for nucleoside identification, the chromatographic (or separation) conditions can now be more easily optimized around assay needs (e.g., speed or elution order) without requiring additional method development steps based on known nucleoside standards.

## Materials and Methods

Cytidine (C), 3-methylcytidine (m<sup>3</sup>C), N<sup>4</sup>-methylcytidine (m<sup>4</sup>C), 5-methylcytidine (m<sup>5</sup>C), adenosine (A), 1methyladenosine (m $^{1}$ A), 2-methyladenosine (m $^{2}$ A),  $N^{6}$ methyladenosine (m<sup>6</sup>A), 8-methyladensine (m<sup>8</sup>A), guanosine (G), 1-methylguanosine (m<sup>1</sup>G), 2-methylguanosine (m<sup>2</sup>G), 7methylguanosine (m<sup>7</sup>G), uridine (U), 3-methyluridine (m<sup>3</sup>U), 5-methyluridine (m<sup>5</sup>U), 2-thiouridine (s<sup>2</sup>U), and 4-thiouridine (s<sup>4</sup>U) were purchased from Carbosynth (Compton, UK). Escherichia coli MRE 600 tRNA, LC-MS grade ammonium acetate, ammonium formate, and ammonium bicarbonate, along with nuclease P1, were purchased from Sigma Aldrich (St. Louis, MO). LC-MS grade water was obtained from Alfa Aesar (Haverhill, MA). LC-MS grade acetonitrile was acquired from Honeywell Burdick & Jackson (Morris Plain, NJ). LC-MS grade formic acid was obtained from Fisher Scientific (Hampton, NH). Phosphodiesterase I and alkaline phosphatase were purchased from Worthington (Lakewood, NJ).

#### LC-MS/MS Analysis of Nucleosides

Analysis of nucleosides was carried out by LC-MS/MS. An ultra-high-performance liquid chromatography (UHPLC) system (Vanguish Flex Quaternary, Thermo Fisher Scientific, San Jose, CA) was used. Reversed phase chromatography was carried out using a high-strength silica column (Acquity UPLC HSS T3, 1.8 µm, 1.0 mm × 50 mm, Waters, Milford, MA). Three different mobile phase conditions were studied. LC method 1 was composed of 5.3 mM ammonium formate in water, pH 4.5, as mobile phase A (MPA), and a mixture of acetonitrile/water (40:60) with 5.3 mM ammonium formate as mobile phase B (MPB). LC method 2 consisted of  $0.1\% v v^{-1}$ formic acid in water, pH 3.0, as MPA, and a mixture of acetonitrile/water (40:60) with  $0.1\% v v^{-1}$  formic acid as MPB. LC method 3 employed 5.3 mM ammonium acetate in water, pH 5.3, as MPA, and a mixture of acetonitrile/water (40:60) with 5.3 mM ammonium acetate as MPB.

The gradient program employed in all three mobile phase conditions consisted of the following: 0% B (from 0 to 6.3 min), 2% B at 13.1 min, 3% B at 16 min, 5% B at 21.4 min, 25% B at 24.6 min, 50% B at 26.9 min, 75% B at 30.2 min (hold for 0.3 min), 99% B at 33 min (hold for 6 min), then returning to 0% B at 39 min. After that, a re-equilibration step at 0% B for 16 min was employed prior to the next injection. A flow rate of  $60~\mu L$  min $^{-1}$  was used. The column temperature was set at  $30~^{\circ}C$ .

Individual nucleosides were resuspended in MPA at a concentration of 20 ng  $\mu L^{-1}$ .  $m^3 C$ ,  $m^4 C$ , and  $m^5 C$  were also prepared in MPA at 5 ng  $\mu L^{-1}$ . Additionally, a mixture of  $m^4 C$  and  $m^5 C$  (20 ng  $\mu L^{-1}$  each) was prepared in MPA. Four micrograms of *E. coli* tRNA was hydrolyzed to nucleosides as previously reported [25]. The solutions were individually analyzed by injecting 1  $\mu L$  (standard nucleosides) or 5  $\mu L$  (*E. coli* hydrolysate) on column.

An Orbitrap Fusion Lumos Tribrid mass spectrometer (Thermo Fisher Scientific) interfaced with a heated electrospray (H-ESI) source (Thermo Fisher Scientific) was used. The analyses were carried out in positive polarity. Full scan data was acquired at a resolution of 120,000, mass range 220–900 m/z, automatic gain control (AGC)  $7.5 \times 10^4$ , and injection time (IT) 100 ms. Data-dependent top speed MS/ MS (1 s cycle) was acquired at a resolution of 15,000, AGC  $1.0 \times 10^4$ , and IT 200 ms. Depending on the experimental goal, CID or HCD was employed for MS/MS fragmentation. The collision energy setting for CID may range from 0 to 100%, while for HCD it may range from 0 to 200 arbitrary units. The other instrumental conditions were the following: quadrupole isolation of 1.5 m/z; ion funnel radiofrequency (RF) level 35%; sheath gas, auxiliary gas, and sweep gas of 30, 10, and 0 arbitrary units, respectively; ion transfer tube temperature of 289 °C; vaporizer temperature of 92 °C; and spray voltage of 3.5 kV.

Data processing was done using the Qual browser of Xcalibur 3.0. Accurate m/z values of the precursor and fragment ions were computed by ChemCalc (http://www.chemcalc.org).

#### Results and Discussion

Model Isomeric Nucleoside Fragmentation in CID and HCD

Reliable and reproducible chromatographic conditions are a requirement for retention time-based differentiation of positional isomers of base-modified nucleosides, due to their identical MH<sup>+</sup> and BH<sub>2</sub><sup>+</sup> ions during mass spectrometry [6]. We initially evaluated CID conditions within the orbitrap to explore the possibility of finding product ions (single or a pattern of multiple) that could uniquely be applied to a specific positional isomer. A change of relative collision energy (CE) from 40 to 100% during CID did not affect the fragmentation pattern as the nucleoside precursor ions yielded nucleobase product ions without exhibiting any other fragment ions (data not shown).

Based on the known features of HCD [18, 19], we reasoned that the HCD capability to trap low mass ions could enable their detection following fragmentation of the nucleobase in multiple pathways at optimal collisional energies. At lower CEs, the HCD-based tandem mass spectrum for ribonucleosides was indistinguishable from that of CID-based tandem mass spectrum (CID-like MS/MS spectra). This observation suggests that HCD can be used for acquiring CID-like MS/MS spectra through conventional N-glycosidic bond cleavage. However, at a CE of 80, low m/z fragment ions were detected, which could potentially be originating from the nucleobase in the HCD-based MS/MS spectrum (e.g., m<sup>5</sup>C in Supplemental Fig. S1). Thus, HCD may offer the potential for both conventional, CID-like MS/MS analysis of modified nucleosides along with additional, higher-energy ring fragmentation that can be used for differentiating structural isomers.

To establish the potential utility of HCD to provide information beyond CID, we initially examined the fragmentation pattern of positional isomers of base-methylated cytidine (m<sup>3</sup>C, m<sup>4</sup>C, and m<sup>5</sup>C) and their reproducible nature in higher-energy HCD-based MS/MS spectra (Table 1). The observed variance in the relative ion abundance (RIA) of product ions was computed for each of positional isomers. The mean and relative standard deviations (RSDs) calculated for the RIAs of nonunique product ions at m/z 126.0667, 109.0402, 108.0562, 83.0609, 82.0293, 81.0453, 56.0500, and 54.0344 are shown in Table 1. These data include replicate analysis (n = 5) at different times of the same day for a fixed LC condition (i.e., LC method 1) and fixed amount of nucleoside injected on column (20 ng). Intraday study data (n = 20) were also compared against data acquired in two different days employing LC method 1 or LC method 2. The effect of sample amount injected on column was also investigated (i.e., 5 or 20 ng).

The examination of the data acquired under these various conditions revealed that the product ions with an RIA  $\geq$  15% exhibited low RSDs ( $\leq$  12%) while the product ions present at lower levels (RIA of < 15%) were less reproducible (RSD  $\leq$  29%). Based on these reproducibility findings, an RIA minimum threshold of 15% was identified as appropriate for development of an HCD-based nucleoside product ion fingerprinting approach. Further, these results demonstrate that this HCD-based approach should be relatively unaffected by experimental parameters such as chromatographic method, analyte retention factor, sample amount, and precursor ion abundance.

Evaluation of the HCD MS/MS spectra of  $\rm m^3C$ ,  $\rm m^4C$ , and  $\rm m^5C$  revealed that most of the observed product ions in the current study correspond well with those reported by Jensen et al. [16] in their CID-based MS³ analysis of the methylated cytosine base ion (BH2+, m/z 126.1). For example, they reported the product ions for  $\rm m^3C$  (at m/z 109.1, 95.1, 83.1, and 69.2),  $\rm m^4C$  (at m/z 109.1, 108.1, 95.1, and 83.2), and  $\rm m^5C$  (at m/z 109.1, 108.1, 83.4, and 81.4). Almost all of these product ions exhibited an RIA of  $\geq$  15% in the HCD MS/MS spectra in this work, except the product ions at m/z 83.1 and 69.2 for  $\rm m^3C$ . These observations suggest the pseudo-MS¹ nature of the HCD spectra of nucleosides.

**Table 1.** RIAs of the Non-Unique Ions with m/z 126.0662, 109.0307, 108.0557, 83.0605, 82.0288, 81.0448, 56.0494, and 54.0337 Observed in the HCD MS/MS Spectra of  $m^3$ C,  $m^4$ C, and  $m^5$ C. Mean and RSD Were Calculated by Accounting Data Acquired Intraday Through the Analysis of Replicates (n = 5), and Interday (2 Different Days) Where Different Amounts of Nucleosides (i.e., 5 and 20 ng) and Different LC Conditions (i.e., LC methods 1 and 2) Were Employed for the Analysis (n = 20). m/z Represented Are Based on Theoretical Calculations

		Intraday ( $n = 5$	5)		Interday and different conditions ( $n = 20$ )				
		m <sup>3</sup> C	m <sup>4</sup> C	m <sup>5</sup> C	m <sup>3</sup> C	m <sup>4</sup> C	m <sup>5</sup> C		
Non-unique ion m/z	126.0667	$100 \pm 0.5$	$100 \pm 0.0$	$100 \pm 0.0$	99 ± 2.3	$100 \pm 0.0$	$100 \pm 0.0$		
	109.0402	$97 \pm 3.5$	$62 \pm 6.2$	$85 \pm 6.2$	$98 \pm 2.6$	$63 \pm 5.8$	$83 \pm 4.6$		
	108.0562	$2.3 \pm 19$	$63 \pm 4.1$	$36 \pm 8.6$	$2.2 \pm 19$	$63 \pm 4.4$	$36 \pm 5.5$		
	83.0609	$6.1 \pm 7.5$	$47 \pm 5.0$	$45 \pm 7.0$	$6.3 \pm 6.9$	$46 \pm 5.8$	$44 \pm 5.6$		
	82.0293	$4.8 \pm 10$	$3.4 \pm 8.2$	$9.5 \pm 15$	$4.7 \pm 16$	$3.3 \pm 13$	$9.9 \pm 11$		
	81.0453	$3.4 \pm 10$	$7.4 \pm 9.7$	$26 \pm 11$	$3.3 \pm 16$	$7.2 \pm 11$	$26 \pm 9.9$		
	56.0500	$2.9 \pm 8.5$	$29 \pm 7.0$	$16 \pm 12$	$2.5 \pm 18$	$29 \pm 8.8$	$16 \pm 10$		
	54.0344	$1.5 \pm 18$	$2.1 \pm 29$	$9.5 \pm 17$	$1.7 \pm 26$	$2.1 \pm 25$	$8.8 \pm 14$		

We next examined the HCD MS/MS data to determine whether product ions are generated that are unique to each isomer. For this purpose, the product ions generated from different positional isomers of methylated cytidine are compared against each other and against the unmodified form (Fig. 1). HCD of m<sup>3</sup>C did not generate any unique product ions at the 15% RIA threshold limit. However, a product ion at m/z 95.0245 is commonly observed for both m<sup>3</sup>C and m<sup>4</sup>C, allowing them to be differentiated from m<sup>5</sup>C. HCD of m<sup>4</sup>C exhibited a unique product ion at m/z 66.0344, allowing it to be differentiated from m<sup>3</sup>C. HCD of m<sup>5</sup>C generated a unique product ion at m/z 81.0453, which distinguishes it from m<sup>3</sup>C and m<sup>4</sup>C. Thus, m<sup>3</sup>C, m<sup>4</sup>C, and m<sup>5</sup>C can be differentiated from each other by HCD based on the presence/absence of specific low mass product ions (Supplemental Table S1). When combined with the other low mass product ions for each isomer discussed above, it appears that HCD can generate isomerspecific product ion fingerprints that can be used for nucleoside identification.

Altering the collisional energy (CE, e.g., 60, 80, 100) during HCD resulted in corresponding alteration of the RIA of observed product ions during MS/MS. However, their abundance is negatively correlated with the abundance of methylated base. In other words, as the energy increases, the abundance of BH<sub>2</sub><sup>+</sup> ion decreased indicating that the product ions are derived from nucleobase. An illustrative comparison of the HCD MS/MS spectra obtained at different CE (60–100, with increments of 20) for m<sup>5</sup>C is shown in Supplemental Fig. S2. Based, in part, on these studies, a CE of 80 was found to provide the best tradeoff between BH<sub>2</sub><sup>+</sup> ion abundance and low mass product ions for ribonucleoside identification.

## HCD Fingerprinting of Model Nucleoside Isomer Mixtures During LC-MS/MS

With the evidence that HCD-based fingerprints could differentiate the methylated cytidine positional isomers, we next examined the applicability of this approach during LC-MS/MS, when such isomers may co-elute. Under typical ammonium acetate- or ammonium formate-based mobile phases, m<sup>4</sup>C and m<sup>5</sup>C are challenging to resolve chromatographically (Fig. 2a).

Conventional CID during LC-MS/MS generates product ion spectra that consist of only the molecular and base ions, which cannot be used to differentiate between these two isomers (Fig. 2b). However, when HCD was used, the fingerprinting approach could be used to identify the presence (or absence) of these isomers (Fig. 2c). The product ion spectra for 100% m<sup>4</sup>C or m<sup>5</sup>C matched those described above. When a mixture of both isomers was analyzed, the product ion spectra reflected the proportion of the individual isomer present in the mixture. In fact, there appears to be a quantitative correlation between the RIAs of the product ions present in the HCD spectra and the relative amount of precursor present in the mixture. Future investigations into using product ion RIAs for reporting initial nucleoside abundance are planned.

# HCD Fingerprinting of Additional Methylated or Other Positional Isomers

To ensure that the HCD fingerprinting approach is generally applicable to modified nucleosides, we examined a few other methylated positional isomers. Methylated adenosines (m<sup>1</sup>A, m<sup>2</sup>A, m<sup>6</sup>A, and m<sup>8</sup>A; Supplemental Fig. S3), methylated guanosines (m<sup>1</sup>G, m<sup>2</sup>G, and m<sup>7</sup>G; Supplemental Fig. S4), and methylated uridines (m<sup>3</sup>U, and m<sup>5</sup>U; Supplemental Fig. S5) were characterized in a manner similar to that described for the methylated cytidines. Table 2 lists the various product ions and RIA for methylated adenosines, guanosines, and uridines. Based on these findings, it appears that this HCD approach is broadly applicable to methylated ribonucleosides.

The differentiation of 2'-O-methylated nucleosides from nucleobase methylated nucleosides is easily accomplished during traditional CID-based or low-energy HCD-based MS/MS by measuring the fragmentation leading to the nucleobase. For the former, a loss of 146 Da (2'-O-methylated ribose) and detection of canonical base is seen, while in the latter the loss of 132 Da (ribose sugar) will yield the methylated base. However, di-methylated nucleosides (e.g., m<sup>4</sup>Cm and m<sup>5</sup>Cm) cannot be differentiated based on CID-type data and would be additional candidates for this HCD-based fingerprinting approach. We tested this possibility by examining the mass spectral fingerprint of m<sup>4</sup>Cm and m<sup>5</sup>Cm using HCD-MS/MS




Figure 1. Fingerprint MS/MS of (a) cytidine, (b) m<sup>3</sup>C, (c) m<sup>4</sup>C, and (d) m<sup>5</sup>C acquired by HCD at CE 80. See Supplemental Table S1 for more information on the peaks labeled

(Fig. 3). The fingerprint of m<sup>4</sup>Cm corresponded with m<sup>4</sup>C, while that of m<sup>5</sup>Cm matched well with m<sup>5</sup>C. This data is consistent with the previous pseudo-MS<sup>n</sup> behavior of HCD demonstrating that HCD can be used to differentiate nucleobase positional isomers, even when the corresponding ribose group is also modified.

Although the data presented above shows the differentiation of methylated isomers by HCD-based MS/MS fingerprinting, we reasoned that positional isomers of other modifications should also be amenable to this approach. To test this possibility, the mass spectral fingerprints of the positional isomers 2- and 4-thiouridine ( $s^2U$  and  $s^4U$ ) were compared (Fig. 4). As anticipated,  $s^4U$  exhibited unique product ions (with RIA  $\geq 15\%$ ) at m/z 129.0122, 86.0064, 68.9799, and 58.9955, and  $s^2U$  yielded unique product ions

at m/z 83.9908, 70.0293 and 59.9908 (Supplemental Table S1). Thus, even for other types of positional isomers, HCD generates unique product ion fingerprints that can be used to precisely identify the parent isomer present in the sample.

## Application of HCD Fingerprinting for LC-MS/MS Analysis of Positional Isomers in Biological Sample

Using a well-defined set of model modified nucleosides, HCD fingerprinting was found to generate reproducible product ions and RIAs, with no observable retention factor effect. To confirm these findings in a more typical laboratory analysis, we applied this approach for the analysis of modified nucleosides obtained from *E. coli* total tRNA. The modified nucleosides in




Figure 2. Extracted ion chromatograms for m/z 258.1090 acquired for (a)  $m^4C$  standard, (b) a 1:1 mixture of  $m^4C$  and  $m^5C$ , and (c)  $m^5C$  standard. Similar chromatographic behavior is observed if LC methods 1 or 3 are used. CID MS/MS of (d)  $m^4C$  standard, (e) 1:1  $m^4C/m^5C$  mixture, and (f)  $m^5C$  standard. HCD-based fingerprint of (g)  $m^4C$  standard, (h) 1:1  $m^4C/m^5C$  mixture, and (i)  $m^5C$  standard

this sample have been well characterized previously [3, 6], thus this sample should provide a good test for the capabilities of HCD-based fingerprinting.

A comparison of the fingerprints acquired for some of the positional isomers detected in this sample is presented in Fig. 5. The observed fingerprint product ion pattern of




Fig. 2. continued.

**Table 2.** Compilation of Peaks m/z and RIAs of Ions Present in the HCD-Based MS/MS Spectra of Positional Isomers of Methylated Adenosines, Guanosines, and Uridines Acquired at CE 80, LC Method 1, and 20 ng of Nucleoside Standard on Column. m/z Represented Are Based on Theoretical Calculations. RIAs Are Represented as Mean (n = 3), with RSD < 10%

$m^1A$	$m^2A$		$m^6A$		m <sup>8</sup> A		$m^1G$		$m^2G$		$m^7G$	$m^3 U$		$m^5U$	
Peak <i>m/z</i> RIA 150.0780 100 133.0514 20 109.0514 20	Peak m/z 150.0780 133.0514 92.0249	RIA 62 100 16	Peak <i>m/z</i> 150.0780 123.0671 94.0405	RIA 100 27 23	Peak <i>m/z</i> 151.0620 150.0780 133.0514 108.0562 106.0405	RIA 51 100 92 17 25	Peak m/z 167.0569 166.0729 153.0412 149.0463 128.0460 121.0514 110.0354 109.0514 67.0296	46 22 16 15	128.0460	RIA 100 21 22 22 15	Peak <i>m/z</i> 167.0569 166.0729 149.0463 124.0511	Peak <i>m/z</i> 127.0508 96.0086	RIA 26 100	Peak m/z 127.0508 110.0242 109.0402 84.0449 82.0293 57.0340 56.0500 54.0340	100




Figure 3. Fingerprint MS/MS of the nucleoside standards (a) m<sup>4</sup>C, (b) m<sup>4</sup>Cm, (c) m<sup>5</sup>C, and (d) m<sup>5</sup>Cm acquired by HCD and CE 80

isomers present in *E. coli* is identical to the specific positional isomer standard's fingerprint described above. As seen in Fig. 5a, the presence of the characteristic peaks at m/z 126.0667 and 109.0402 in the MS/MS spectrum of the

precursor with m/z 258.1090 indicated the presence of a methylcytidine modification. Additionally, the absence of peaks at m/z 95.0245 and 66.0344, and the presence of the unique ion at m/z 81.0453, confirmed the identity of the




Figure 4. Fingerprint MS/MS of (a) uridine, (b) s<sup>2</sup>U, and (c) s<sup>4</sup>U acquired by HCD at CE 80. See Supplemental Table S1 for more information on the peaks labeled




Figure 5. Fingerprint MS/MS of (a) m<sup>5</sup>C and (b) m<sup>4</sup>Cm detected in commercial E. coli tRNA

m<sup>5</sup>C positional isomer. In addition to the characteristic unique ions, the RIAs of non-unique and unique product ions matched well with those specified in Table 1 and the standard m<sup>5</sup>C shown in Fig. 1d. The same can be concluded

for other positional isomers detected in *E. coli* tRNA (Supplemental Fig. S6). Thus, the fingerprints acquired from biological samples can be correlated to the fingerprint acquired for specific standards through spectral matching.




Figure 6. Comparison of extracted ion chromatograms (XIC) and CID MS/MS at CE 40 and HCD MS/MS at CE 80, of Am, m<sup>2</sup>A, and m<sup>6</sup>A detected in *E. coli* under different LC conditions. (a) XIC for *m/z* 282.1202 employing LC method 3. (b) XIC for *m/z* 282.1202 employing LC method 2. CID MS/MS spectra of m<sup>2</sup>A for LC methods 3 and 2 are shown in (c) and (g), respectively. CID MS/MS spectra for m<sup>6</sup>A for LC methods 3 and 2 are shown in (e) and (i), respectively. HCD MS/MS spectra for m<sup>6</sup>A for LC methods 3 and 2 are shown in (f) and (j), respectively.

The fingerprint for m<sup>4</sup>Cm detected during this analysis is shown in Fig. 5b. For this 2'-O-modification, its presence could be confirmed by the occurrence of characteristic ions and RIAs, identical to the respective non-ribose-modified methylcytidine (i.e., m<sup>4</sup>C). Detection of this nucleoside (as well as m<sup>5</sup>C) in the commercial E. coli tRNA sample was unexpected as neither of these modifications have previously been reported in E. coli tRNAs. These findings suggest that either these modifications are present in tRNA locations not previously reported or, more likely, the sample might be contaminated with (degraded) ribosomal RNA (rRNA), as both modifications have been previously reported in E. coli rRNA [26, 27] and electrophoresis of the commercial tRNA sample did not show any intact rRNA (data not shown). Nevertheless, the data shows that even low abundant modifications can be detected and differentiated by HCD-based fingerprinting during standard bioanalytical characterization of RNA nucleoside digests.

In addition to identifying positional isomers when a standard LC method is ineffective at separating those compounds, another advantage of using HCD fingerprinting is when the LC method is modified (e.g., for method improvement or to enhance MS performance). To illustrate the retention factor insensitivity of HCD fingerprinting, the same E. coli total tRNA nucleoside digest analyzed above was also analyzed using an alternative LC method [9] consisting of  $0.1\% v v^{-1}$  formic acid in water, pH 3.0, as MPA, and acetonitrile/water 40:60,  $0.1\% v v^{-1}$  formic acid, as MPB (LC method 2). The extracted ion chromatograms (XICs) presented in Fig. 6 show three of the nucleosides detected at m/z 282.1202 (i.e., Am, m<sup>2</sup>A, and m<sup>6</sup>A). Using an ammonium acetate mobile phase, the three peaks, from left to right, correspond to Am, m<sup>2</sup>A, and m<sup>6</sup>A, respectively (Fig. 6a). However, using formic acid in the mobile phase, the elution order for m<sup>2</sup>A and Am shifts (Fig. 6b). HCD fingerprinting confirms that it is m<sup>2</sup>A, rather than m<sup>6</sup>A, which shifts under these mobile phase conditions, and the fingerprint for both positional isomers is essentially identical to those obtained using the ammonium acetate-based mobile phase (Fig. 6c-j). Thus, based on these HCD fingerprints, one can readily identify nucleobase isomers without resorting to additional analyses using standards.




Figure 7. Comparison of extracted ion chromatograms (XICs) for  $\Psi$ ,  $t^6A$ , oQ and  $k^2C$  detected in *E. coli* using (a) LC method 1 and (b) LC method 2. \* indicates  $t^6A$  diastereomer. Corresponding HCD fingerprints from each of the different LC methods for  $\Psi$  (c) method 1/(d) method 2,  $t^6A$  (e) method 1/(f) method 2, oQ (g) method 1/(h) method 2, and  $t^2C$  (i) method 1/(i) method 2




Fig. 7. continued.

## HCD-Based Fingerprinting for Any Modified Ribonucleoside

While an important utility of HCD is to use product ion fingerprints to differentiate positional isomers, the more general case holds that any nucleoside presents a characteristic and reproducible fingerprint. For instance, when an E. coli hydrolysate was analyzed by employing either LC method 1 or 2, retention time differences varied differently for each modified nucleoside (Fig. 7a, b). Nevertheless, the HCD-based fingerprints of the ribonucleosides were highly reproducible (Fig. 7c–j). This observation is illustrated for the RNA modifications pseudouridine ( $\Psi$ ), epoxyqueousine (oQ),  $N^6$ -threonylcarbamoyladenosine ( $t^6$ A), and 2-lysidine ( $t^2$ C).

As employed, HCD-based fingerprints cannot be used to differentiate diastereomers such as D-allo- and L-t<sup>6</sup>A. The latter is a naturally occurring RNA modification (known simply as t<sup>6</sup>A) present in tRNA of a variety of organisms (bacteria, fungi, plants, and protists) [28]. The former has been recently shown by Matuszewski et al. [29] to be present in *E. coli* tRNA and is linked to the hydrolysis of cyclic t<sup>6</sup>A (ct<sup>6</sup>A) under the mild basic conditions employed for RNA hydrolysis. Both diastereomers can be observed in the XICs shown in Fig. 7a, b, in which D-allo-t<sup>6</sup>A (indicated by an asterisk) is an earlier eluter when compared to L-t<sup>6</sup>A. This new approach, however, could be used to provide additional validation or strong support

for diastereomers when nucleosides have different elution times but share all mass spectral features, including the HCD fingerprint, in common.

#### **Conclusions**

We have found that HCD-based dissociation of modified ribonucleosides during standard LC-MS/MS analysis generates nucleoside-specific product ion fingerprints. These fingerprints can be used in a standard constant-neutral loss mode technique for nucleoside identification. Importantly, positional isomers of modified nucleosides are readily differentiated based on these fingerprints. Moreover, even in those cases where positional isomers co-elute, individual fingerprints can be deconvoluted from the mixture data to identify which isomers are present. These fingerprints are insensitive to chromatographic conditions, precursor ion abundance, and sample matrix.

This HCD-based fingerprinting approach enables the detection and confirmation of RNA modifications based exclusively on their mass spectrometric behavior, which allows the possibility of developing an automated protocol for the detection of RNA modifications though spectral matching. Automated nucleoside detection would be possible via the compilation of a spectral library of HCD fingerprints acquired for nucleoside standards and/or nucleosides detected in different organisms,

which would then facilitate matching against analytes detected in biological samples. Additionally, this work reveals that the pseudo-MS<sup>n</sup> nature of HCD generates a rich collection of purine and pyrimidine ring product ions, which should prove useful in situations where a modified nucleoside is not previously characterized.

## **Funding Information**

Financial support of this work was provided by the NSF (CHE1507357), DTRA (HDTRA1-15-1-0033), NIH (GM061115 to P.A.B. and OD018485 to P.A.L.) and an NIH T32 Training Grant (GM113770) supporting Cody M. Palumbo. The generous support of the Rieveschl Eminent Scholar Endowment and the University of Cincinnati for these studies is also appreciated.

## References

- Boccaletto, P., Machnicka, M.A., Purta, E., Piątkowski, P., Bagiński, B., Wirecki, T.K., de Crécy-Lagard, V., Ross, R., Limbach, P.A., Kotter, A., Helm, M., Bujnicki, J.M.: MODOMICS: a database of RNA modification pathways. 2017 update. Nucleic Acids Res. 46, D303–D307 (2017)
- Miyauchi, K., Kimura, S., Suzuki, T.: A cyclic form of N6threonylcarbamoyladenosine as a widely distributed tRNA hypermodification. Nat. Chem. Biol. 9, 105–111 (2012)
- Kellner, S., Neumann, J., Rosenkranz, D., Lebedeva, S., Ketting, R.F., Zischler, H., Schneider, D., Helm, M.: Profiling of RNA modifications by multiplexed stable isotope labelling. Chem. Commun. 50, 3516–3518 (2014)
- Kang, B.-i., Miyauchi, K., Matuszewski, M., D'Almeida, G.S., Rubio, M.A.T., Alfonzo, J.D., Inoue, K., Sakaguchi, Y., Suzuki, T., Sochacka, E., Suzuki, T.: Identification of 2-methylthio cyclic N6threonylcarbamoyladenosine (ms2ct6A) as a novel RNA modification at position 37 of tRNAs. Nucleic Acids Res. 45, 2124–2136 (2017)
- Zallot, R., Ross, R., Chen, W.-H., Bruner, S.D., Limbach, P.A., de Crécy-Lagard, V.: Identification of a novel epoxyqueuosine reductase family by comparative genomics. ACS Chem. Biol. 12, 844–851 (2017)
- Russell, S.P., Limbach, P.A.: Evaluating the reproducibility of quantifying modified nucleosides from ribonucleic acids by LC-UV-MS. J. Chromatogr. B. 923-924, 74-82 (2013)
- Su, D., Chan, C.T.Y., Gu, C., Lim, K.S., Chionh, Y.H., McBee, M.E., Russell, B.S., Babu, I.R., Begley, T.J., Dedon, P.C.: Quantitative analysis of ribonucleoside modifications in tRNA by HPLC-coupled mass spectrometry. Nat. Protoc. 9, 828–841 (2014)
- Thüring, K., Schmid, K., Keller, P., Helm, M.: Analysis of RNA modifications by liquid chromatography—tandem mass spectrometry. Methods. 107, 48–56 (2016)
- Basanta-Sanchez, M., Temple, S., Ansari, S.A., D'Amico, A., Agris, P.F.:
  Attomole quantification and global profile of RNA modifications:
  epitranscriptome of human neural stem cells. Nucleic Acids Res. 44,
  e26-e26 (2016)
- Heiss, M., Reichle, V.F., Kellner, S.: Observing the fate of tRNA and its modifications by nucleic acid isotope labeling mass spectrometry: NAIL-MS. RNA Biol. 14, 1260–1268 (2017)
- Pomerantz, S.C., McCloskey, J.A.: Analysis of RNA hydrolyzates by liquid chromatography-mass spectrometry. Methods Enzymol. 193, 796– 824 (1990)
- Chan, C.T.Y., Pang, Y.L.J., Deng, W., Babu, I.R., Dyavaiah, M., Begley, T.J., Dedon, P.C.: Reprogramming of tRNA modifications controls the

- oxidative stress response by codon-biased translation of proteins. Nat. Commun. **3**(937), (2012)
- Cai, W.M., Chionh, Y.H., Hia, F., Gu, C., Kellner, S., McBee, M.E., Ng, C.S., Pang, Y.L.J., Prestwich, E.G., Lim, K.S., Ramesh Babu, I., Begley, T.J., Dedon, P.C.: A platform for discovery and quantification of modified ribonucleosides in RNA: application to stressinduced reprogramming of tRNA modifications. Methods Enzymol. 560, 29–71 (2015)
- Sakaguchi, Y., Miyauchi, K., Kang, B.-i., Suzuki, T.: Nucleoside analysis by hydrophilic interaction liquid chromatography coupled with mass spectrometry. Methods Enzymol. 560, 19–28 (2015)
- Squires, J.E., Preiss, T.: Function and detection of 5-methylcytosine in eukaryotic RNA. Epigenomics. 2, 709–715 (2010)
- Jensen, S.S., Ariza, X., Nielsen, P., Vilarrasa, J., Kirpekar, F.: Collisioninduced dissociation of cytidine and its derivatives. J. Mass Spectrom. 42, 49–57 (2007)
- Rose, R.E., Quinn, R., Sayre, J.L., Fabris, D.: Profiling ribonucleotide modifications at full-transcriptome level: a step toward MS-based epitranscriptomics. RNA. 21, 1361–1374 (2015)
- Olsen, J.V., Macek, B., Lange, O., Makarov, A., Homing, S., Mann, M.: Higher-energy C-trap dissociation for peptide modification analysis. Nat. Meth. 4, 709–712 (2007)
- Senko, M.W., Remes, P.M., Canterbury, J.D., Mathur, R., Song, Q., Eliuk, S.M., Mullen, C., Earley, L., Hardman, M., Blethrow, J.D., Bui, H., Specht, A., Lange, O., Denisov, E., Makarov, A., Horning, S., Zabrouskov, V.: Novel parallelized quadrupole/linear ion trap/orbitrap tribrid mass spectrometer improving proteome coverage and peptide identification rates. Anal. Chem. 85, 11710–11714 (2013)
- de Graaf, E.L., Altelaar, A.F.M., van Breukelen, B., Mohammed, S., Heck, A.J.R.: Improving SRM assay development: a global comparison between triple quadrupole, ion trap, and higher energy CID peptide fragmentation spectra. J. Proteome Res. 10, 4334–4341 (2011)
- Abdelhameed, A.S., Kadi, A.A., Abdel-Aziz, H.A., Angawi, R.F., Attwa, M.W., Al-Rashood, K.A.: Multistage fragmentation of ion trap mass spectrometry system and pseudo-MS3 of triple quadrupole mass spectrometry characterize certain (E)-3-(dimethylamino)-1-arylprop-2-en-1ones: a comparative study. Sci. World J. 2014, 1–14 (2014)
- McCloskey, J.A., Graham, D.E., Zhou, S., Crain, P.F., Ibba, M., Konisky, J., Söll, D., Olsen, G.J.: Post-transcriptional modification in archaeal tRNAs: identities and phylogenetic relations of nucleotides from mesophilic and hyperthermophilic Methanococcales. Nucleic Acids Res. 29, 4699–4706 (2001)
- Xu, L., Liu, X., Sheng, N., Oo, K.S., Liang, J., Chionh, Y.H., Xu, J., Ye, F., Gao, Y.-G., Dedon, P.C., Fu, X.-Y.: Three distinct 3-methylcytidine (m3C) methyltransferases modify tRNA and mRNA in mice and humans. J. Biol. Chem. 292, 14695–14703 (2017)
- Asano, K., Suzuki, T., Saito, A., Wei, F.-Y., Ikeuchi, Y., Numata, T., Tanaka, R., Yamane, Y., Yamamoto, T., Goto, T., Kishita, Y., Murayama, K., Ohtake, A., Okazaki, Y., Tomizawa, K., Sakaguchi, Y., Suzuki, T.: Metabolic and chemical regulation of tRNA modification associated with taurine deficiency and human disease. Nucleic Acids Res. 46, 1565–1583 (2018)
- Ross, R., Cao, X., Yu, N., Limbach, P.A.: Sequence mapping of transfer RNA chemical modifications by liquid chromatography tandem mass spectrometry. Methods. 107, 73–78 (2016)
- Kimura, S., Suzuki, T.: Fine-tuning of the ribosomal decoding center by conserved methyl-modifications in the Escherichia coli 16S rRNA. Nucleic Acids Res. 38, 1341–1352 (2010)
- Andersen, N.M., Douthwaite, S.: YebU is a m5C methyltransferase specific for 16 S rRNA nucleotide 1407. J. Mol. Biol. 359, 777–786 (2006)
- Yacoubi, B.E., Bailly, M., Crécy-Lagard, V.d.: Biosynthesis and function of posttranscriptional modifications of transfer RNAs. Ann. Rev. Genet. 46, 69–95 (2012)
- Matuszewski, M., Wojciechowski, J., Miyauchi, K., Gdaniec, Z., Wolf, W.M., Suzuki, T., Sochacka, E.: A hydantoin isoform of cyclic N6threonylcarbamoyladenosine (ct6A) is present in tRNAs. Nucleic Acids Res. 45, 2137–2149 (2017)